



tecno-lab

ATTREZZATURE SCIENTIFICHE
DA LABORATORIO & ASSISTENZA TECNICA
Strumentazione Consumabili



tecno-lab s.r.l.
di Galesi e Pistoni

Via L. Abbiati, 22/A-B • 25131 Brescia - Italia
Tel. 030.3582505 r.a. • Fax 030.3582517
www.tecnolab.bs.it • info@tecnolab.bs.it

KS 1000 Scrubber

Fumes Neutralization

Designed to neutralize toxic fumes produced during mineralization and safeguard lab operators without requiring any connection to tap water



Space saving and ease-of-use

The KS 1000 Scrubber guarantees maximum workplace safety by neutralizing corrosive and toxic fumes produced during Kjeldahl digestion in two steps: dilution and neutralization.

- High-efficiency gas scrubbing for maximum safety in the laboratory
- Tap water-free for reduced water consumption and no need for a dedicated water supply
- A small footprint for a high space saving
- Designed for DK and DKL connection
- 4 adjustable suction power levels
- Intuitive and easy-to-use
- High-quality construction materials
- Usable for both acid and basic neutralization



INSTRUMENT - CODE

KS 1000 Scrubber	100-240 V / 50-60 Hz	F307A0660
------------------	----------------------	-----------

TECHNICAL DATA

CONSTRUCTION MATERIAL	Metal structure and PVC components
WORKING IN CONTINUOUS	Admitted
FLOW REGULATION	Up to 25 L/min
VACUUM REGULATION	Down to 200 mbar abs
CAPACITY OF GLASS TANKS	3.5 L
POWER INPUT	130 W
CURRENT CONSUMPTION	500 mA
WEIGHT	11 Kg 24.3 lb
DIMENSIONS (WXHxD)	315x335x360 mm 12.4x13.2x14.2 in