

MULTIVAP 6 CONCENTRATOR BASIC SOLUTION

The **MultiVap 6** is a basic, cost effective sample preparation instrument featuring 6 independent channels. Without transferring, samples are directly concentrated in the vial thereby reducing the risk of transfer loss, improving recovery and saving precious time.

The nitrogen flow rate of each channel can be adjusted independently as well as the needles that can be adjusted up and down.

The aluminum alloy heating block guarantees temperature uniformity for all samples.



TECHNICAL FEATURES

- Possibility to run up to 6 samples of 40ml
- Parallel concentration model
- Possibility to simultaneously or separately work with six channels
- Nitrogen blow in combination with heating concentration mode
- Dry block heating technology to avoid potential water vapor interference
- P.I.D. temperature control with digital display
- Possibility to choose the best position of nitrogen needles to ensure the most efficient evaporation
- Dedicated nitrogen flow regulator for each sample channel
- Compatibility with collection vials of Sepaths UP
- Unique “open view” to control the sample volume level during the evaporation process
- Temperature range: ambient to 120°C
- Weight: 8kg
- Dimensions: (LxWxH cm) 13x35x20
- Power supply voltage: 230V 50/60Hz

Included

- N° 6 - 40ml Collection bottles
- N° 1 - 220V Power Cord
- N° 1 Two-way connector O.D4mm/O.D6mm
- N° 1 PU Tubing O.D4mm, 2m
- N° 1 PU Tubing O.D6mm, 2m
- N° 6 Nitrogen needles
- N° 1 Flange tool

Ordering information U-SP9005-C01

ACCESSORIES

P/N	DESCRIPTION	P/N	DESCRIPTION
LVS40C-ET	40ml sample collection vial - flat bottom	AWZ7P0	Aluminium block Tube sizes 2ml (12mm OD) & 16x100mm
Y02060017	40ml Graduated collection vial - 1ml end point (1pz)		
B4VRS0	Aluminium block - Tube sizes 15mm & 22mm OD	AWZ7Q0	Aluminium block Tube sizes 2ml (12mm OD) & 16x100mm
AWZ7M0	Aluminium block - Tube size 18mm OD		

■ *Specially designed adapters available*