



tecno-lab

ATTREZZATURE SCIENTIFICHE
DA LABORATORIO E ASSISTENZA TECNICA
Strumentazione Consumabili



tecno-lab s.r.l.
di Galesi e Pistoni

Via L. Abbiati, 22/A-B • 25131 Brescia - Italia
Tel. 030.3582505 r.a. • Fax 030.3582517
www.tecnolab.bs.it • info@tecnolab.bs.it

RC 2

Heating Plates

Two-position heating plates with aluminum top



Practical and Easy to use

Heating plates with aluminum top and temperature regulation. Designed for the longest lifespan and to guarantee excellent resistance to chemicals. It accepts different combinations of AluBlocks™ that can maximize temperature control and save handling time.

- Maximum temperature 370 °C
- Accepts flasks of up to 15 L volume
- Uniform heat distribution with coated aluminum alloy heating plate



INSTRUMENT - CODE

RC2	230 V / 50-60 Hz	F20700430
RC2	115 V / 50-60 Hz	F20710430

TECHNICAL DATA

	RC2
HOUSING MATERIAL	Epoxy painted aluminum structure
PLATE MATERIAL	Aluminum alloy coated with a special protection
PLATE DIMENSION	135 mm
MODELS	2-position
TEMPERATURE RANGE	Room temp. - 370 °C
POWER INPUT	1200 W
WEIGHT	3.3 Kg (7.3 lb)
DIMENSIONS (WxHxD)	340x90x246 mm (13.4x3.5x9.7 in)