



tecno-lab

ATTREZZATURE SCIENTIFICHE
DA LABORATORIO E ASSISTENZA TECNICA
Strumentazione Consumabili



tecno-lab s.r.l.
di Galesi e Pistoni

Via L. Abbiati, 22/A-B • 25131 Brescia - Italia
Tel. 030.3582505 r.a. • Fax 030.3582517
www.tecnolab.bs.it • info@tecnolab.bs.it

RC 5

Heating Plates

Single-position heating plate with aluminum top



Excellent Heating Performance

The RC 5 is equipped with a 135 mm diameter aluminum alloy top plate ensuring a great temperature homogeneity and optimum heat transfer across the entire surface.

Designed to guarantee ease of use, the longest lifespan and excellent resistance to chemicals.

- Maximum temperature 310 °C
- Accepts flasks of up to 15 L volume
- Hot plate warning indicators
- IP 42 certified



INSTRUMENT - CODE

RC 5	230 V / 50-60 Hz	F20700610
RC 5	115 V / 50-60 Hz	F20710610

TECHNICAL DATA

	RC 5
HOUSING MATERIAL	Epoxy painted aluminum structure
PLATE MATERIAL	Aluminum alloy
PLATE DIMENSION	135 mm
TEMPERATURE RANGE	Room temp. - 310 °C
HOT TOP WARNING SYSTEM	Illuminated Icon (when the temperature is above 50 °C)
POWER INPUT	630 W
PERMISSIBLE AMBIENT TEMPERATURE	5 - 40 °C
PROTECTION CLASS	IP 42
WEIGHT	1.4 Kg (3.1 lb)
DIMENSIONS (WxHxD)	160x85x270 mm (6,3x3,3x10,6 in)